

Digital Medical Scope VersaCamTM





Digital Solution for Viewing, Showing, and Recording Various Clinical Images

The VersaCam™ provides the digital documentation of clinical images that a physician examines.



Wide Range of Interchangeable Lenses

Interchangeable lens adapters allows the use in various medical applications such as an ophthalmoscope, anterior eye segment scope, otoscope, and dermatoscope.








Improved Informed Consent



Captured images can be displayed directly on the 3.5-inch color LCD monitor with a 2x magnification or on the sub monitor*, which facilitates easier patient education and informed consent.

* AV output of VersaCam $^{\!\mathsf{TM}}$ conforms with NTSC format.







Real time display on sub monitor

Easy and Quick Image Filing

Captured images can be imported to image filing software. As one of which, NAVIS-EX is available. About 80,000 images can be saved on a micro SD memory card.







Versatile Digital Medical Scope



High Adaptability

This handheld medical scope achieves ease of use for image capture of a patient who is seated or in the supine position. This unit is also attachable to a slit lamp.







Excellent Portability

The VersaCam™ is easily portable, weighing only 400 g including the ophthalmoscope adapter. A lithium-ion battery provides up to 3 hours of continuous operation.



Easy One-handed Operation

The user friendly design enables the easy one-handed operation.




The still image mode and video mode can be switched.

An image can be captured by pressing the button or foot switch.



High Quality Digital Imaging

The clear-color camera with 40° angle of view achieves high quality digital imaging.



The focus can be adjusted with the dial to quickly capture images at various ranges even on the retina or tympanic membrane.

Various Lens Adapters





Clear-color digital images with 40° angle of view can be captured and saved.



Eye cup

Peripheral Devices (optional)








Foot switch

The slit lamp adapter is useful to attach the unit to the slit lamp, which allows easier operation. Using the adapter with the foot switch minimizes blurring of images due to hand movement.


Clinical Images



Normal eye




Geographic atrophy



Neonatal eye

Anterior Eye Segment Scope

With this adapter, the flash color can be switched from white to blue which is useful to evaluate fluorescein staining.










This adapter allows imaging of the ear canal and tympanic membrane with a wide angle of view.





Polarization is a property of light that oscillates in all directions. Filtering the light reflected from surface with the polarized light (PL) filter is helpful for inner skin evaluation.

An image of pigmented lesion such as skin disorders and moles can be captured with this adapter. With the PL filter, pigment distribution on the inner skin can be observed.





Tympanic Membra

Without PL filter

With PL filter

General Lens







This adapter can be operated easily similar to a general camera and used for clinical images of the throat, skin and other tissues. The LED flash is helpful for capturing images in dimly lit environments.

Endoscope Adaptor

With this adapter, an endoscope*1*2 can be attached to the VersaCam™.

- *1 NIDEK does not offer endoscopes.
- *2 Please contact a NIDEK distributor to confirm if the endoscopes can be attached to this adapter.



VersaCam™ DS-10 Specifications

M	ain	body

Main body	
Camera Resolution	1,920 x 1,080 pixels
File format	JPEG (still image), H.264 (video)
Memory	8GB micro SD card
Display	3.5-inch color LCD
Interface	Mini USB, AV output (phone jack, NTSC)
Power supply	Lithium-ion rechargeable battery
	AC 100 to 240 V
	50 / 60 Hz
Dimensions / Mass	89 (W) x 44 (D) x 205 (H) mm / 254 g
	3.5 (W) x 1.7 (D) x 8.1 (H)" / 0.56 lbs.

Charging station

Interface	Mini USB
Dimensions / Mass	49 (W) x 156 (D) x 83 (H) mm / 333 g
	1.9 (W) x 6.1 (D) x 3.3 (H)" / 0.73 lbs.

Lens adapter (optional)

Ophthalmoscope	
Туре	Non-mydriatic
Angle of view	40°
Dioptric compensation	-20 to +20 D
for patient's eyes	
Light source	Infrared LED (observation),
	White LED (observation and capturing)
Dimensions / Mass	135 (L) x 44 (Dia.) mm / 140 g
	5.3 (L) x 1.7 (Dia.)" / 0.31 lbs.
Anterior Eye Segment Scope	
Angle of view	41° (diagonal)
Light source	White LED, Blue LED
Dimensions / Mass	45.2 (L) x 44.5 (Dia.) mm / 85 g
	1.8 (L) x 1.7 (Dia.)" / 0.19 lbs.
Otoscope	
Focusing range	5 to 50 mm
Light source	White LED
Dimensions / Mass	69 (L) x 44 (Dia.) mm / 68 g
	2.7 (L) x 1.7 (Dia.)" / 0.15 lbs.
Dermatoscope	
Image capture range	10 mm (diagonal)
Optical modes	Polarization, Non-polarization
Light source	White LED
Dimensions / Mass	56 (L) x 44 (Dia.) mm / 113 g
	2.2 (L) x 1.7 (Dia.)" / 0.25 lbs.
General Lens	
Angle of view	88° (diagonal)
Light source	White LED
Dimensions / Mass	25 (L) x 44 (Dia.) mm / 49 g
	1.0 (L) x 1.7 (Dia.)" / 0.11 lbs.
Endoscope adapter	
Dimensions / Mass	50 (W) × 18 (D) × 59 (H) mm / 34 g
	2.0 (W) x 0.7 (D) x 2.3 (H)" / 0.07 lbs.

	Interface	Phone Jack
	Dimensions / Mass	73 (W) x 102 (D) x 31 (H) mm / 265 g
		2.9 (W) x 4.0 (D) x 1.2 (H)" / 0.58 lbs.



Slit lamp adapter (optional)

	3.0 (W) x 5.2 (D) x min. 8.0 (H)" / 1.9 lbs.
Dimensions / Mass	// (W) x 132 (D) x min. 202 (H) mm / 850 g



Manufacturer: MedimagingIntegrated Solution Inc.

1F, No. 7, R&D Rd II, Hsinchu Science Park, Hsinchu, Taiwan 30076 EC Rep: Obelis S.A.: Bd. General Wahis 53, 1030, Brussels, Belgium

Product / Model name: Digital Medical Scope DS-10EA

Digital Medical Scope DS-10F Digital Medical Scope DS-10A Digital Medical Scope DS-10E Digital Medical Scope DS-10D Digital Medical Scope DS-10G

Caution: U.S. Federal Law restricts this device to sale, distribution, and use by or on the order of a physician or other licensed eye care practitioner.

Specifications may vary depending on circumstances in each country.

Specifications and design are subject to change without notice.



HEAD OFFICE (International Div.) 34-14 Maehama, Hiroishi Gamagori, Aichi 443-0038,

Japan TEL: +81-533-67-8895 URL: http://www.nidek.co.jp [Manufacturer]

TOKYO OFFICE (International Div.)

3F Sumitomo Fudosan Hongo Bldg., 3-22-5 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan TEL: +81-3-5844-2641 URL: http://www.nidek.com

NIDEK INC.

47651 Westinghouse Drive Fremont, CA 94539, U.S.A. TEL: +1-510-226-5700 +1-800-223-9044 (US only) URL: http://usa.nidek.com

NIDEK S.A.

Europarc 13, rue Auguste Perret 94042 Créteil, France TEL: +33-1-49 80 97 97 URL: http://www.nidek.fr **NIDEK TECHNOLOGIES Srl**

Via dell'Artigianato, 6 / A 35020 Albignasego (Padova), Italy TEL: +39 049 8629200 / 8626399 URL: http://www.nidektechnologies.it NIDEK (SHANGHAI) CO., LTD. NIDEK SINGAPORE PTE. LTD. Rm 915, China Venturetech Plaza, No.819 Nanjing West Rd,

Jing An District, Shanghai China 200041 TEL: +86 021-5212-7942 URL: http://www.nidek-china.cn 51 Changi Business Park Central 2 #06-14

The Signature Singapore 486066 TEL: +65 6588 0389

