



### Features

- ▶ Easy to use one button system from 24hrs up to 48hrs of BP measurement
- ▶ The most accurate Hyper, Hypo-tension detection algorithm
- ▶ Useful diagnosis for titrating medication
- ▶ Non-dipping diagnosis for Left Ventricle Hypertrophy
- ▶ Durable & patented cuff with inner pocket for user friendly wearing system  
With adult, large adult, Extra adult cuff sizes
- ▶ Easy data management interface with configuration, analysis, reporting & diagnosis



# LAB24™

24-HR ABP

## 24-48Hour Ambulatory Blood Pressure Monitoring System


### Clinical application

#### White-coat hypertension


- Resistant hypertension
- Masked hypertension
- Childhood hypertension
- Efficacy of anti-hypertensive drug therapy on a 24-hour basis
- Nocturnal hypertension
- Episodic hypertension and/ or anxiety disorders
- Hypotensive symptoms
- Changes in diet and daily routine designed to reduce hypertension



### Easy interpretation Window based SW with USB cable retrieve connection



Hourly average report



ABP date graph with Treshhold tab

### General Specifications

Validations:	BHS (A/A), ESH, ANSI/AAMI SP10:2002
BP Technique:	Oscillometry with step deflation
BP Range:	25-260 mmHg
HR Range:	40-200 bpm
Sampling Periods:	3 adjustable (awake, asleep, and special)
Time Intervals:	5, 10, 15, 20, 30, 45, 60, 90 and 120 min.
Event Key:	Program to allow additional patient readings
Data Storage:	Flash memory stores up to 250 readings Two
Power Source:	AA batteries
Weight:	284 grams, including batteries
Dimensions:	12 x 7 x 3cm
PC Requirements:	USB port, Windows® XP or tater



North America (U.S.A) : Bionet America, Inc.  
Tel : +1-714-734-1760

E-mail : sales@BionetUS.com  
www.BionetUS.com

South Korea : Bionet co., Ltd  
Tel : +82-70-7585-6496(Int'l)  
+82-2-6300-6483(Domestic)

E-mail : sales@ebionet.com(Int'l)  
domestic.sales@ebionet.com(Domestic)  
www.eBionet.com